KENTDENTAL

- fr Instruments endodontiques manuels
- en Endodontic Manual Instruments
- it Strumenti endodontici manuali
- es Instrumentos manuales de endodoncia
- pt Instrumentos endodônticos manuais
- de Manuelle endodontische Instrumente
- nl Handmatige endodontische instrumenten
- bg Ръчни ендодонтски инструменти
- pl Ręczne instrumenty endodontyczne
- ro Instrumente endodontice manuale
- ru Ручные эндодонтические инструменты
- ar أدوات المعالجة اللبية اليدوية لطب الاسنان
- fa ابزارهای دستی درون دندانی

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Instruments endodontiques manuels Limes K, Broches, Racleurs, Tire-nerfs

Description

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K-Files, K-Reamers, Hedstroem Files, Barbed broaches

Description

K files, K-Reamers, Hedstroem Files and barbed broaches are endodontic manual instruments with taper of 1% for barbed broaches and 2% for the other instruments.

The instruments are composed of operative part made of 304 stainless steel (also containing nickel-chromium), a handle made of plastic (polybutylene terephthalate) and a stopper made of silicone. The devices are for single use, supplied non-sterile and must be sterilised by a healthcare professional before use.

Intended Use

K files, K-Reamers Hedstroem Files, and barbed broaches are endodontic instruments, intended to explore, shape, and clean the root canal system during no surgical endodontic treatment.

Users

Endodontic instruments can only be used in a clinical or hospital setting, in accordance with good dental practices, by qualified dental professionals such as general practitioners and endodontists.

Target population

These endodontic instruments are intended for use on adults and children with a tooth affected by pulpitis or pulp necrosis.

Indications

K files, K-Reamers Hedstroem Files, and Barbed broaches are indicated for the non-surgical endodontic treatment.

Contraindications

The product contains Nickel and Chromium and should not be used on patients with known allergic sensitivity to these metals.

Complications

In cases of complex canal anatomy, per-operative risks (fracture, perforation, abutment or other) could occur and lead to a risk of infectious processes.

Characteristics

See table in the appendix.

Precautions for use

- Do not use the instruments when storage conditions are suspected.
- Inspect the packaging before use and do not use the instruments if the packaging is damaged.
- If you have any doubts about the product identity, do not use it.
- These products are sharp, remember to always handle them by the handle, to avoid the risk of injury when handling and using.
- For your own safety, please wear personal protective equipment (gloves, glasses, mask).
- Remember to check the instrument state before use. If the instrument is damaged or shows signs of wear, do not use it.
- Do not force the instrument during use.
- If the instrument does not easily progress, inspect and clean the cutting flutes, then irrigate.

- Check the instrument and the working part frequently during use, inspecting for signs of distortion, elongation or wear, such as uneven flutes, dull spots. These signs indicate that the devices are not able to fulfil the intended use with the required safety level, instruments must be discarded.
- The instrument should not be completely immersed in sodium hypochlorite (NaOCl) solution. Only the active part of the instrument, in contact with the patient, may be immersed in a NaOCl solution with a concentration not exceeding 5% for a maximum of 5 minutes.
- Instrument supplied non-sterile: Perform the cleaning and sterilisation procedure before use, in accordance with the instructions for use.

Warnings

- These products should only be used by properly trained dentists-surgeons.
- These products are limited to dental treatment purposes and are prohibited for uses other than root canal treatment or retreatment.
- Any serious incident involving the device shall be reported to the manufacturer and to the competent authority of the member state in which the user and/or patient is established.
- Strictly adhere to the cleaning - disinfection, sterilisation, and use of the product, to minimize as much as possible, the following risks to patients and/or users:
 - Instrument fracture.
 - Cross-contamination.
 - Ingestion or inhalation of the instrument.
 - Toxic or allergic reactions caused by device residues.
 - Paresthesia due to over instrumentation.

Technical performance

- Respect for root canal anatomy (respect for original canal axis).
- Adapted geometry in the apical zone (taper $\geq 2\%$ and $\geq 1\%$ only for barbed broaches).
- Shaping ability (debris removal, cutting efficiency).

Indirect clinical benefits

- The use of these endodontic instruments allows dentists to preserve the tooth from potential avulsion, by extracting the neurovascular bundle from the root canal(s) and then proceed with their shaping for a filling.
- Treatment with endodontic instruments helps to relieve patients in cases of inflammation of the pulp or other pain or discomfort associated with the progression of caries disease or infection.

Product Storage Conditions

Store the product in a dry and clean place, away from sunlight and moisture.

Steps of use

Preparation before use

The product is not supplied sterile and should be cleaned, disinfected and sterilised before use (thorough cleaning and disinfection are essential for effective sterilisation) by following the instructions for use below.

General recommendations	<ul style="list-style-type: none"> • For all metal devices, it is recommended to use anticorrosion cleaning and disinfecting agents. • For your own safety, wear protective equipment (gloves, safety glasses, mask and gown). • Do not use cleaning or disinfecting agents containing phenol or aldehyde. • Perform all stages involving detergent solutions in a well-ventilated area. Avoid inhalation of the solution and contact with skin and eyes. • Always follow the instructions for use provided by the product manufacturer.
-------------------------	---

Limitations on use	<ul style="list-style-type: none"> The expiration time of the device is 5 years, respect the expiration date of the product indicated on the packaging, do not use the device after this date.
Materials needed	<ul style="list-style-type: none"> Gloves, masks, gown according to the manufacturer's recommendations of the cleaning and disinfecting agent. Running or demineralized water. Detergent Neodisher Mediclean Forte. Small soft brushes. Container. Ultrasonic bath or washer disinfectant. Class B sterilisation apparatus. <p>Remark: All materials used should be cleaned and replaced regularly. Identify materials used for each step of the process (initial treatment, cleaning or rinsing).</p>

1	<p>Initial treatment</p> <p>Immediately place products in a container with tap water (< 40°C) and detergent for 5 to 15 min.</p> <p>Rinse the products with tap water at 20-40°C for 1min.</p> <p>If there is a pause before the next step, make sure the instruments stay moist by packing them in wet wipe. Do not exceed 1 hour waiting time.</p> <p>Remarks:</p> <ul style="list-style-type: none"> Do not use fixing agents or hot water (> 40°C), as this causes fixation of residues. Follow instructions and observe concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or other defects on the instruments).
2	<p>Preparation before cleaning</p> <p>If the devices have visible impurities, it is recommended to manually pre-clean by brushing them under tap water at 20-40°C at least 1 min with a soft brush until the impurities disappear.</p> <p>Remark:</p> <p>Follow instructions and observe concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or other defects on the devices).</p>
3	<p>Visual inspection</p> <p>Inspect the products and discard any damaged products (broken, untwisted, or abnormally bent).</p>

Automatic	4	Cleaning/Rinsing/Drying Use a washer-disinfector and place the instruments in the machine on a tray. Perform a cleaning cycle with detergent. Do not use high-concentration cleaning agents as this may damage the surface of the instrument.	
	5	Automated thermal disinfection in a washer-disinfector: A 5-minute disinfection cycle at 93°C.	
	6	Automated drying: Drying of instruments by the washer/disinfector drying cycle. If necessary, additional manual drying can be carried out using a lint-free towel or sterile compressed air.	
		Remarks: <ul style="list-style-type: none"> Follow the instructions and concentrations given by the manufacturer of the detergent solution. Follow the instructions of washer-disinfector and check that the success criteria after each cycle have been met as stated by the manufacturer. Use only washers/disinfectors that comply with the EN ISO 15883 standard in force and are maintained and validated regularly. 	
	7	Visual inspection After cleaning and disinfecting, make sure the product is visibly clean. Check for any marks on the product (dents, cracks, deformations, scratches or corrosion). Repeat steps 4-5-6 if the product is not visually clean. Discard all damaged products.	
	8	Packaging Place the instruments in a paper-plastic pouch for steam sterilisation in compliance in effect with ISO 11607 and EN 868 standards. Remark: Seal the pouches according to the recommendations of the pouch manufacturer. If a heat-sealer is used, the process must be validated, and the heat-sealer must be calibrated and qualified.	
	9	Sterilisation In accordance with the mandatory French regulations applicable in France (Cycle Prion): Apparatus: Class B. Temperature: 134°C Cycle time: 18 min. Absolute pressure: 2.2 bar. Control the physico-chemical indicators and cycle parameters.	In accordance with the current ISO 17665-1 international standard: Apparatus: Class B. Temperature: 134°C Cycle time: 5 min. Absolute pressure: 2.2 bar.
		Remarks: <ul style="list-style-type: none"> When sterilising multiple instruments in one autoclave cycle, ensure that the steriliser's maximum load is not exceeded. Place the pouches in the steam steriliser according to the recommendations given by the steriliser manufacturer. Use only a pre-vacuum air removal steam steriliser that meets the requirements in effect of EN 13060 (class B, small sterilizer) and EN 285 (standard sterilizer) with saturated steam. Rapid sterilisation or the use of sterilisation methods for unpackaged instruments is not permitted. Do not use a high-pressure steam sterilizer with a temperature exceeding 200 degrees Celsius, including the drying process. 	

10	<p>Storage after sterilisation After sterilisation: Remove the product from the autoclave. Let it cool to room temperature for at least 30 minutes. Do not use additional cooling. Check that the sterilisation wraps or pouches are not damaged. In case of damage, a complete restoration should be carried out. Store the products in a dry, clean and dust-free environment, at temperatures specified for paper-plastic pouches by the steam sterilisation manufacturer.</p> <p>Remark: Refer to the label and instructions for use.</p>
----	--

Protocol of use

- Preliminary confirm the working length based on a radiological judgement of the image.
- Use a rubber barrier system (dam), limiting the risk of ingestion or inhalation of the instruments.
- Remove the pulp with the barbed broaches
- Use a suitable sterile type K-File for catheterization to advance with a slight push of 2-3 mm, reaching the working length one or more times.
- Irrigate carefully with sodium hypochlorite in a concentration not exceeding 5%.
- Gradual use of other types of files with increasingly larger diameters depending on the size of the canal, while always respecting the anatomy of the canal and the working length identified.
- Irrigate carefully with sodium hypochlorite between each file re-entry into the canal.
- When using the instrument, it should be pushed as gently as possible.
- The violent use of instruments is strictly prohibited.
- The use of viscous chelating agents can relatively reduce the release of metal ions from the instrument.

Mode of action

K-files and K-Reamers:

They are used in insertion movement with a clockwise/counter-clockwise rotation of one quarter or one eighth of a turn and then a removal.

Hedstroem Files:

They are used in insertion and removal movements without any rotational action.

Barbed broaches:

They are used in insertion, with a rotation and then a traction

Disposal of the device

After use, to avoid any risk of infection and microbial contamination, the devices must be disposed of in a safety container identical to the one required for sharp, cutting, and disposable piercing instruments, in accordance with good dental practices.

Identifying Symbols

	Manufacturer		Importer
	Date of manufacture		Swiss Representative
	Expiration date		Caution
	European compliance		Unique Device Identifier
	Store away from direct sunlight		Keep dry
	Assortment		Do not use if damaged packaging
	See instructions for use		K Reamers
	Quantity		Stainless steel material
	Non-sterile		Stop made of silicone
	K Files		Endodontic treatment
	H Files		Batch Code
	Catalogue number		Autoclave at the specified temperature
	Do not reuse		Medical device
	Barbed Broaches		

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Dental Addict SUISSE S.A.