38.000

2.0000 | 3.0000 | 3.0000 | 3.0000 | 3.0000 | 3.0000 | 3.0000 | 3.0000 | 3.0000 | 3.0000 | 3.0000 | 3.0000 | 3.0000

FOLD

EN ENGLISH

Product Description
Silventum™ Clinpro™ Clear Fluoride Gel contains a 2.1% sodium fluoride, water-based coating that is applied topically to enamel and dentin for the treatment of hypersensitive teeth.Clinpro Clear fluoride treatment contains 9,000 ppm fluoride as well as added calcium, phosphate and xylitol. This product is supplied in single use, unit-dose LPep packages containing 0.5 ml (0.5 grams) of fluoride coating.

Intended Use
Fluoridated tooth coating.

Indications for Use
Clinpro Clear fluoride treatment is indicated for the treatment of hypersensitive teeth.

Contraindications
None known.

Intended Users
Educated dental professionals who have theoretical and practical knowledge on usage of dental products.

Intended Patient Population
Intended patient population includes all patients requiring dental treatment unless the patient's condition limits the use. The safety and performance of Clinpro™ Clear Fluoride Treatment have not been clinically evaluated in pregnant and lactating women.

Clinical benefit
Treatment of tooth hypersensitivity.

Precautions and Warnings
For Patients
Discontinue use of other prescriptive fluoride preparations (i.e. gels, rinses) for 24 hours following Clinpro Clear fluoride treatment application.

For Dental Personnel
Avoid use in patients that may have known hypersensitivity to the product or its ingredients. To limit ingestion, use only enough product to form a thin coating on the desired treatment area.
Safety Data Sheet (SDS)/Safety Information Sheet (SIS) can be obtained from www.Solventum.com or contact your local subsidiary.

Potential Side Effects
Gastrointestinal upset, local irritations, or local or systemic allergic reactions.

Residual Risks
In situations where the Instructions for Use was not followed or the product was misused (i.e., excess ingestion or aspiration of product, a tendency to swallow exposures the following risks may occur: headache lightheadedness, dizziness, presyncope, toxicity, eye or airway irritation.

Instructions for Use –
General application instructions:
1. For best results, teeth should be clean and free of excess saliva prior to application.
2. Press the single-use LPep package by gently pushing the thumb and index finger to completely squeeze the product out of the reservoir and into the attached dispensing well.
3. Remove the applicator brush from the package by gently pulling the handle.
4. Using the applicator brush, apply the product to the teeth. Reload the brush as needed. It is not necessary to use all the product provided. Use enough product to form a coating on the desired treatment area.
5. After application, do not rinse. You may see a clear coating on the teeth.
6. The fluoride coating will gradually wear away. To achieve the maximum benefit, please instruct the patient to avoid food, beverages, oral rinses for at least 10 minutes after application.

Storage Conditions
This product is designed to be stored and used at room temperature. Shelf life at room temperature is 24 months. Ambient temperatures routinely higher than 27°C/80°F, or lower than 10°C/50°F, may reduce shelf life. Do not freeze or expose to extreme heat. See outer package for expiration date.

Handling after Use
Dispose of the single use only. Properly dispose of used or contaminated product to avoid health risks due to improper handling.

SDS and Disposal
See the Safety Data Sheet/Safety Information Sheet (available at www.Solventum.com or contact your local subsidiary) for disposal and transportation information.

To be disposed in accordance with the local or regional regulations.

Please pay special attention to the disposal of contaminated waste to avoid health risks due to improper handling.

Customer Information
Please report a serious incident that has occurred in relation to the device to Solventum and the local competent authority (EU) or local regulatory authority.

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

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Symbol Glossary

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Authorized Representative in the European Community / European Union	EC REP	Indicates the authorized representative in the European Community / European Union. Source: ISO 15223, 5.1.2; 2014/35/EU; and/or 2014/30/EU
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code	LOT	Indicates the manufacturer's batch code so that the batch or lot can be identified. Source: ISO 15223, 5.1.5
Catalogue number	REF	Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed. Source: ISO 15223, 5.3.7
Do not re-use		Indicates a medical device that is intended for one single use only. Source: ISO 15223, 5.4.2
Consult instructions for use or consult electronic instructions for use		Indicates the need for the user to consult the instructions for use. Source: ISO 15223, 5.4.3
CE Mark	CE	Conformity to all applicable European Union Medical Device Regulations and Directives.
Unit dose L-Pep		Indicates unit-dose L-Pep

For symbols glossary, see eIFU.Solventum.com

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© Solventum 2026. Solventum, le logo S et Clinpro sont des marques commerciales de Solventum ou de ses filiales.