ENGLISH

Intended Use

A dental composite resin for anterior and posterior restorations.

Product Description

Solventum™ Filtek™ Supreme Flowable Restorative is a low viscosity, visible-light activated, radiopaque, flowable nanocomposite. The restorative is packaged in multiple syringes and capules. It is available in a variety of tooth-colored shades. The shades offered with Filtek Supreme Flowable Restorative were designed to coordinate with shades offered with Filtek™ Supreme Ultra, Filtek™ Ultimate, and Filtek™ Supreme XTE Universal Restorative. The resin of this restorative contains bisGMA, TEGDMA and Procrylat resins. The fillers are a non-agglomerated/non-aggregated surface modified 20 nm silica filler, a non-agglomerated/non-aggregated surface modified 75 nm silica filler, and a surface modified zirconia/silica cluster filler (comprised of 20 nm silica and 4 to 11 nm zirconia particles). The aggregate has an average cluster particle size of 0.6 to 10 microns. The inorganic filler loading is approximately 63% by weight (46% by volume).

Solventum™ Filtek™ Supreme Flowable Restorative is a type 1, low viscosity, visible-light activated, radiopaque flowable, complies with ISO 4049:2019 Dentistry – Polymer-based restorative Material and ISO 6874: Class 2 Light Cure Pit and Fissure Sealant.

Capsules are radiopaque, with a value of 1.9 mm of aluminum. Aluminum has a radio-opacity equivalent to that of 1.9 mm of material having a radio-opacity equivalent to 1 mm of aluminum and 2 mm of aluminum is opaque to x-rays.

Indications for Use

- Direct anterior and posterior restorations (including occlusal surfaces)
- Base/liner under the restorative directes
- Repair of metal-ceramic and ceramic-based direct and indirect restorations
- Pit and fissure sealant

Intended Patient Population

Patients requiring dental treatment unless the patient's medical condition is a contraindication.

Contraindications

None known

Warnings and Precautions

Educated dental professionals who have theoretical and practical knowledge on the usage of dental products.

Clinical benefit

Restoration of oral esthetics and function.

Precautions and Warnings

For Dental Professional

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce this risk, avoid contact with uncured material. Use gloves to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves during placement and finishing is recommended. Acrylates may penetrate commonly used gloves. If product contacts glove, remove and discard glove, wash hands immediately with soap and water. If uncured material enters the eye, medical attention is needed.

Risks and Undesirable Effects

- Allergic reaction (i.e., redness, swelling, itching, burning sensation near the site of application)
- Post-operative sensitivity/pain

Interactions

Prescription Use

1. **Prophy:** Teeth should be cleaned with pumice and water to remove surface stains.
2. **Shade Selection:** Prior to isolation of tooth, determine desired shade(s) of Filtek Supreme Flowable Restorative using a standard VITAPAN® classical shade guide. Shade selection accuracy may be enhanced via a mock-up; place some of the material adjacent to the prepared tooth. Moisture will change the appearance of the shade, so evaluate thickness and site of the restoration. Cure Evaluate the shade match under different lighting sources. Remove the mock-up material with a plastic instrument. Do not use a sharp instrument as this may damage the tooth surface. If achieved, for sealants, a shade guide is not required. Select a shade that is close in color.
3. **Isolation:** A rubber dam is the preferred method of isolation. Cotton rolls and an evacuator can also be used.

Directions

1. **Matrix band:** Formation of tight proximal contacts is accomplished solely through careful matrix technique. Wedging is needed to provide a tight gingival seal and to ensure a separation between teeth to compensate for thickness of the matrix band. Where possible, use of a predontoured sectional matrix band is preferred. However, a thin circumferential matrix band may be used. Filtek™ Supreme Flowable Restorative is not to be used to establish proximal contacts. Filtek Supreme Flowable Restorative should not hold the proximal matrix band in place.
2. **Pulp protection:** If a pulp exposure has occurred and the situation warrants a direct pulp capping procedure, use a minimum amount of calcium hydroxide liner. Follow the manufacturer's instructions for the application of Vitrebond™ Plus Light Cure Glass Ionomer. Vitrebond Plus liner/bases may also be used to line areas of deep cavity excavation. Filtek Supreme Flowable Restorative should not be used with eugenol containing materials, as the cure may be inhibited.
3. **Adhesive System:** To bond Filtek Supreme Flowable Restorative to tooth structure, use of a Solventum™ dental adhesive system is recommended. Refer to adhesive system instructions for complete directions and precautions for the products.

 After curing the adhesive, continue to maintain isolation from blood, saliva and other fluids and proceed immediately to placement of Filtek Supreme Flowable Restorative.
4. **Cure:** Filtek Supreme Flowable Restorative is a light-cure composite. Refer to adhesive application.
5. **Warming:** Syringes can be warmed up to 37°C for one minute in commercially available composite warmer devices. Capsules can be warmed up to 70°C for one hour in commercially available composite warmer devices.
6. **Dispensing:**

 Syringe dispensing: Filtek Supreme Flowable Restorative can be delivered directly from the dispensing tip. Protective eyewear for patients and staff is recommended when using the dispensing tip.

 5.1. Remove the syringe tip cap. Attach a new dispensing tip. The dispensing tip need to be discarded and replaced with each new patient.

 5.2. Twist the disposable dispensing tip securely onto the syringe. As needed, bend the metal cannula in any direction, up to a 90° angle to facilitate access to the site.

 5.3. Holding the tip away from the patient and any dental staff, express a small amount of Filtek Supreme Flowable Restorative to ensure that the material flows freely from the tip.

 5.4. If plugged, remove the dispensing tip and clean the syringe opening. Replace dispensing tip and again express composite. Filtek Supreme Flowable Restorative is provided unit a dispensing gun equipped with an appropriate instrument.

 5.5. Before first use of the syringe during treatment, a standard dental protective sleeve should be used. Place the functional syringe and the protective sleeve in a syringe holder. Push the plunger of the syringe down to the end, then remove the plunger, opening the cannula for use. Using a protective sleeve helps to avoid contamination of the syringe between patients. See the syringe holder instructions for further details.

 Single-Dose Capsule: Insert capsule into Filtek™ Restoratives Dispenser. Refer to separate restorative dispenser instructions for full instructions and precautions on dispensing and delivery of composite.
7. **Placement:**

 7.1. Avoid intense light in the working field. Exposure of paste to intense light may cause premature polymerization.

 7.2. Dispense the desired amount of material. Avoid dispensing more than is needed for the preparation, holding the tip of the syringe away from the preparation. Dispensing more than is needed may cause the tip to adhere to the material. For proximal areas, hold the tip against the matrix to aid material flow into the proximal box.

 7.3. Place and light cure restorative in increments of 2 mm or less. Refer to step 4 for curing.
8. **Curing:** This product is intended to be cured by exposure to a halogen or LED light with a minimum intensity of 550 mW/cm² in the 400-500 nm range. Use of lower intensity curing lights or curing lights with a limited light source, such as a Solventum™ curing light. Hold the light guide tip as close to the restorative as possible during light exposure.

Cure Time		
Shades	**Increment Depth**	**All halogen light's and LED lights (with output 1000 – 2000 mW/cm²)**
Opaque Shades	1,5 mm	40 sec.
	2,0 mm	20 sec.
All Other Shades		10 sec.

8. **Complete the restoration:**

 8.1. Finishing Line application: Place a composite restorative material, such as a Filtek™ restorative material directly over the cured Filtek Supreme Flowable Restorative. Follow the manufacturer's instructions regarding placement, curing, finishing, occlusal adjustment and polishing.

 8.2. Pit and Fissure Sealants: Gently remove the inhibited layer remaining after light curing with slurry of pumice or non-fluoride prophy paste.

 8.3. Direct Restorative Applications:

 8.3.1. Contour restoration surfaces with fine finishing diamonds, burs or stones.

 8.3.2. Use a composite with an appropriate shade and texture to replace the enamel layer. Refer to manufacturer's instructions for placement, curing, finishing, occlusal adjust occlusion by removing material with a fine polishing diamond or stone.

 8.3.3. Polish restoration with composite polishing instrument such as Sof-Lex™ finishing system.

Storage Conditions

Do not store unsafely inspect the product for damage. If in doubt, please contact Solventum or your local distributor

This product is designed to be used at room temperature. If stored cold, the syringe can be warmed in a commercial composite warmer prior to use (at no higher than 70°C / 158°F, for no longer than 1 hour), up to 25 times. The capsule may be warmed in a commercial composite warmer prior to use (at no higher than 70°C / 158°F, for no longer than 1 hour).

Store in a cool, dry location allow product to reach room temperature prior to use. Keep away from direct light. 18 months. Ambient temperatures above 25°C (77°F) may reduce shelf life. See outer package for expiration date.

Do not expose restorative materials to elevated temperatures, or to intense light.

Do not store materials in proximity to eugenol containing products.

After use of dispensing syringe, remove and discard syringe tip. Clean syringe body as needed. **Discard used syringe.tip and.sleeve.in.appropriate waste.stream.. Do not reuse syringe sleeve or tip.**

Handling after Use

Use new, unused disposable gloves during all clinical applications for the syringe to reduce the risk of contamination.

Place the disposable tip on the syringe and place the protective sleeve over the syringe. Use the tip to punch through the end of the protective sleeve.

Avoid contact between the reusable parts (e.g., the body of the syringe) and the patient's mouth.

Use only a light intensity that will allow product to reach room temperature prior to use. Keep away from direct light.

18 months. Ambient temperatures above 25°C (77°F) may reduce shelf life. See outer package for expiration date.

Do not expose restorative materials to elevated temperatures, or to intense light.

Do not store materials in proximity to eugenol containing products.

After use of dispensing syringe, remove and discard syringe tip. Clean syringe body as needed. **Discard used syringe.tip and.sleeve.in.appropriate waste.stream..**

Handling after Use

Use new, unused disposable gloves during all clinical applications for the syringe to reduce the risk of contamination.

Place the disposable tip on the syringe and place the protective sleeve over the syringe. Use the tip to punch through the end of the protective sleeve.

Avoid contact between the reusable parts (e.g., the body of the syringe) and the patient's mouth.

Use only a light intensity that will allow product to reach room temperature prior to use. Keep away from direct light.

Disposal of the syringe unless it is sterilized and remains in the unopened and undamaged blister pack.

Dispose of in accordance with local regulations for dental offices and/or hospitals during the use and reprocessing of the device.

The capsules are for single use only. Properly dispose of used or contaminated capsules to avoid health risks because of improper handling.

Please refer to the Filtek™ Restoratives Dispenser IFU for cleaning, disinfection and sterilization instructions.

If the syringe is not visibly clean, i.e. if there is accidental contamination of the syringe with blood or saliva, discard the syringe immediately. Do not place the syringe into a disinfection bath or a cleaning/disinfection device (washer-disinfector).

Ingredients

Silane Treated Ceramic, Silane Treated Silica, Zirconia/Silica Cluster, Trimethylene Glycol Dimethacrylate, Ytterbium Fluoride, Reacted Polycaprolactone Polymer, Diphenyldiindonium Hexafluorophosphate

SDS/SS

- Safety Data Sheet (available at www.Solventum.com)
- See SDS/SS for disposal information or to be disposed in accordance with the local or regional regulations.
- Please pay special attention to the disposal of contaminated waste to avoid health risks due to improper handling.

Customer Information

- Please report a serious incident occurring in relation to the device to Solventum and the local competent authority (EU) or local regulatory authority.
- Summary of the Safety and Clinical Performance can be found at https://ec.europa.eu/tools/eudamed. Basic UDI-DI for capsules: 07979988401020000000115NN; basic UDI-DI for syringes: 07979988401020000000115NL

Warranty

Solventum warrants this product will be free from defects in material and manufacture at the time of purchase. **EXCEPT AS STATED IN THE PRECEDING SENTENCE, SOLVENTUM MAKES NO OTHER WARRANTY OR GUARANTY OF FITNESS FOR A PARTICULAR PURPOSE, MERCHANTABILITY OR OTHERWISE COMMERCIAL IN CHARACTER.** If this product is defective within the warranty period, your exclusive remedy and Solventum's sole obligation shall be repair or replacement of the Solventum product.

Limitation of Liability

Except where prohibited by law, Solventum will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For more information, see eIFU.Solventum.com

Symbol Title	Symbol	Symbol Glossary Description and Reference
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Authorized Representative in European Community / European Union	EC REP	Indicates the authorized representative in the European Community / European Union. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code	LOT	Indicates the manufacturer's batch code so that the batch or lot can be identified. Source: ISO 15223, 5.1.5
Catalogue number	REF	Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed. Source: ISO 15223, 5.3.7
Do not re-use		Indicates a medical device that is intended for one single use only. Source: ISO 15223, 5.4.2
Medical Device	MD	Indicates the item is a medical device. Source: ISO 15223, 5.7.7
Unique device identifier	UDI	Indicates a carrier that contains unique device identifier information. Source: ISO 15223, 5.7.10
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8
Swiss Authorised Representative	CH REP	Indicates the authorised representative in Switzerland. Source: Swissmedic.ch
CE Mark	CE 0123	Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
Refer to instruction manual/ booklet		To signify that the instruction manual/booklet must be read. ISO 7010 M002

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℗ Solventum 2025. Solventum, the S logo, Filtek, Vitrebond, and Sof-Lex are trademarks of Solventum or its affiliates.
Solventum US LLC
Solventum Germany GmbH